DR. RAJENDRA PRASAD CENTRE FOR OPHTHALMIC SCIENCES All India Institute of Medical Sciences Ansari Nagar, New Delhi-29

Ref. No. 03/SSO(RPC)/RS/PAC/2019-20

Subject: Purchase of Surgical Microscope with Intra-Operative Spectral Domain OCT - 01 No. for Dr. R.P.Centre, AIIMS, New Delhi-29 on proprietary basis- <u>Inviting</u> comments thereon.

proposal submitted by M/s Carl Zeiss (India) Bangalore Pvt. Ltd., (Authorized representative of As per decision taken/ approved by Competent Authority of Dr. R.P.Centre AIIMS for the purchase of subject cited equipment from M/s. Carl Ziess, Germany on proprietary basis. The Carl Ziess, Germany) and PAC certifications are attached & uploaded on website.

Dr. within 15 days from the date of issue/uploading of the notification giving reference R.P.Centre at AlIMS on or before 07.10.2019 upto 12.30 P.M., failing which it will be comments, if any, from any manufacturer regarding proprietary nature of the equipment/item presumed that any other vendor is having no comment to offer and case will be decided on The above documents are being uploaded for open information to submit objections/ 03/SSO(RPC)/RS/PAC/2019-20. The comments should be sent to Sr. Stores Officer, merits.

Encl: Related documents enclosed. 1. PAC Certificate enclosed.

SR. STORES OFFICER (RPC)

To Consider purchase of Intraoperative OCT integrated operating microscope with HD Camera. Compatible with ngenuity 3D Visual on proprietary basis

Main Microscope:

Aprochomatic optics with anti-reflex multi coating

Page | 4

- Motorizled zoom system with zoom ratio 1:6 magnification factors: 0.4x-2.4x
 - Focussing range 70 mm
- Speed control for zoom and focus
- Tiltable binocular tube f= 170 mm, with Integrated image inverter.
- Pair of high eyepoint widefield eye pieces 10 xwith diopter setting from -8D to +5D
- Apochromatic objective f= 200 mm with carrier ring.

 Total magnifications: 4.3x to 25.5x with eyepiece 12.5x and objective lens f=200mm Field of view: 8.6 mm to 51.8 mm with eyepiece 12.5x and objective lens f=200mm
 - Integrated Slit illumination; Slit width 0.2, 2.0, 3.0, 4.0mm & SlIt height 12mm
 - beam splitter should be integrated in the microscope body
- HD camera should be integrated in the microscope body without any external attachment. Camera controls unit should be integrated in the stand.

Built-in assistant's Microscope:

- magnification to achieve magnification for main surgeon & assistant, independent fine focusing Integrated Assistant microscope with electrical zoom magnification, with programmable
- Inclinable Binocular tube with integrated image inverter.
- SCI (Stereo coaxial illumination) for constant brilliance and brightness, red reflex illumination and surrounding field illumination both are adjustable.
 - Pair of high eye point wide field eye pieces 10x with diopter setting from -8D to +5D,
 - Provision of red reflex for assistant with equal brightness

Coupling

- Range of adjustment 60 mm x 60 mm. Control of automatic reset of XY movements.
- Provision of inversion of XY direction of travel via foot control, Speed control for XY.

Illumination

- SCI (Stereo coaxial illumination) for constant brilliance and brightness, red reflex illumination and surrounding field illumination both are adjustable.
 - Fiber light guide, Integrated Xenon illumination system with 180W xenon lamp with back up lamp 180W xenon with availability of Halogen filtered illumination.
 - Integrated 408 nm UV filter for protection against infrared exposure
 - Blue Blocking Filter, Provision of retina protection device
- Provision of system of magnetic clutches for all locks for positioning of microscope across surgical field

Intraoperative OCT:

- Spectral Domain OCT wavelength 840nm
- Scanning speed 27000 A scan per second
 - Axil resolution 5.5micron.



Scan modes for Live: 1 Line, 5 Line & Cross hair. Scan modes for capture: 1 Line, 5 line & Cube. Touch screen control through Callisto eye.

OCT Module should be integrated in Microscope body without any external attachment.

OCT Module & Microscope must be from same manufacturer.

All the components of OCT should be factory integrated in microscope body & Stand.

Wide Angle viewing system:

Non contact, Autoclavable wide angle viewing system.

Non contact lenses 60D & 128D - 2 sets.

Markerlesstoric IOL alignment

Markerlesstoric IOL alignment, Assistance Markerless License* Reference Axis, Rehixs, Main Incision & Paracenthesis, LRI, Z ALIGN® Marker based & Markerless

Target axis for toric IOL alignment
Can be set directly on CALLISTO eye 3.6 or imported via from IOLMaster 700*
One or three lines, Position relative to yellow reference axis
OCT scans 2.9mm for High resolution OCT scan & 5.8mm for large overview to visualize & access graft orientation.

Distortion free computer enhanced intraoperative OCT image to visualize detailed structures in the correct physiological shape & better view/ observe the Irido-corneal angle.

Floor Stand

Magnetic clutches for effortless movement and positioning, Built in maneuvering handles

Facility to change to back up lamp in event of lamp failure by fast action change Lamp intensity adjustment via foot control panel

Progressive speed adjustments

Wireless programmable 14 function foot control panel.

Storage facility of magnification, motor speed, configuration of foot control panel, lamp brightness and focal plane for at least 9 different users

Facility for non sterile release of suspension arm

Accessories

HD Video Recorder should be integrated in microscope stand.

IDIS facility- Superimpose of OCT image in eye pieces.

Should be compatible with Calisto Imaging System
** It should be US FDA/ European CE approved.

ALL INDIA INSTITUTE OF MEDICAL SCIENCES ANSARI NAGAR, NEW DELHI 110029

PROPERIETARY/SPECIFIC BRAND GOODS CERTIFICATE

Item/Type/Model No. Required along with specification.

Microscop Intraoperative Operating Demain OCT Surgical with

Spectra

or accessory for an existing equipment. Is the item a spare part attachment ri

: No (New Equipment)

Name of the manufacturers/suppliers of the item proposed by the indenter. ei

: M/s. Carl Ziess, Germany

Are they sole manufacturers/ sole distributors of the item. 4

: Sole Manufacturer

the same can't be procured. Demanding equivalent specification available in the recommended item from those offered officer should bring out comparative advantages/ cost effectiveness of the market to meet the job requirement envisaged. If the answer is yes, why Is there any other item with similar/ from others. vi

S.

What were the efforts made to locate alternative source of supply or use 6

the not

Why open/limited tender can't be other substitutes. 1

As per the best of our knowledge other specifications are with available required manufacturer.

> resorted to, for locating alternative source.

The required specification are not other the with manufacturer available . .

> rates are reasonable or not. Are the proprietary items 00

Yes, prices are reasonable

Any other justification for procuring item from single source. 6

Integrated HD Video Recorder is not available with the other facility of Inbuilt OCT, manufacturer The

I certify that the item at Sr. No. 1 above is required to be procured on single tender basis as the source of supply is definitely known/ the specific brand proposed was advantageous in meeting our functional requirements and limited tender system could be dispensed with as they would serve no useful purpose in this particular case.

(Strike out whichever is not applicable)

Signature of Indenter (Demanding officer)

(Head of Department)



Carl Zeiss Meditec AG 10589 Berlin

Dr. R. P. Centre For Ophthalmic Sciences All India Institute of Medical Sciences New Delhi - 110029 Ansari Nagar, The Chief,

Carl Zeiss Meditec AG
Berlin Ste
Max-Dohrn-Strasse 8-10
10589 Berlin, Germany

Phone: +49 (0) 30/ 854001-301

e-mail: juoka.ruhanen@zeiss.com +49 (0) 30/ 854001-123 Fax

PROPRIETARY CERTIFICATE FOR ZEISS OPMI LUMERA 700 SURGICAL MICROSCOPE AND RESCAN 700 INTEGRATED INTRAOPERATIVE OPTICAL COHERENCE TOMOGRAPHY (OCT) SYSTEM

Dear Sir.

We hereby certify that the OPMI LUMERA 700 and RESCAN 700 AIO, the first surgical microscope ophthalmic surgery with integrated intraoperative optical coherence tomography (OCT) system, is the proprietary product of Carl Zeiss Meditec AG, Germany Equipped with ZEISS RESCAN 700, ZEISS OPMI LUMERA 700 takes surgical microscopy to a whole new level with integrated intraoperative OCT. The ophthalmic surgeons visualize transparent structures of the anterior and posterior segments directly in the eyepieces of the surgical microscope; see exactly where they are scanning with the scan location marker and move the scan independently of the surgical microscope. ZEISS RESCAN 700 gives the ophthalmic surgeons more information during retinal or corneal surgery in the eyepieces of the surgical microscope, allowing them to see structures in new ways; helping them back up their decisions, thereby improving surgical techniques or outcomes without compromising their surgical workflow.

Address of Record: Goeschwitzer Strasse 51 - 52 07745, Jena, Germany

Address for Deimery.
Car Zeiss Meditec AG
Max-Dohm-Stracse 8 - 10 10589 Berlin, Germany

Banks: Deutsche Gank Jena Account: 6245,36900 (BIC 870 700 00) JBAN DE90 8707 0000 0624 5369 00 BIC SWIFT: DEUT DE 8EXXX

Account. 258072809 (\$1C 820 409 00) IRAN: DE-11 8204 0000 11758 0728 00 BICL SWIFT, COBADCIFYXX

Commercial Register: Local Court Jens HIRB 205623

VAT-4D No., DE 811 922 737 WEEF-Reg -No.; DE55298748

Charman of the Supe Dr. Michael Kaschke

Board of Management.
Dr. Ludwin Monz (CEO)
Autus Felia Wehmer
Am Willem de Ger



Carl Zeiss Meditec AG Page 2 of 2

System has best-in-class and unique features such as:

CALLISTO eye with Markerless toric IOL alignment
 ASSISTANT funcions for Rhexis, Main Incision & Paracenthesis, LRI, Z ALJGN (toric IOL)
 Reference image/axis from IOLMaster
 Guided image quality check for optimizing light intensity, magnification & centration prior to matching the live image to the reference axis image

- OCT scans with 2.9mm depth for high resolution images as well as OCT scans with 5.8mm depth for images providing an excellent and large overview (e.g. to visualize corneal graft orientation).
- Undistorted visualization (correct uspect ratio) of computer enhanced intraoperative OCT images of the eye m
- OCT scans in posterior segment for vitreo-retinal surgery, are not only supported with RESIGHT but also in combination with a flat contact glass or indirect contact glass.

Best regards

Carl Zeiss Meditec AG

Anonevery

Head of Marketing and Business Management Ana Neves

Jukka Ruhanen Senior Business Manager APAC